

For North East Hampshire and Farnham Place and Surrey Heath Place only

## North East Hampshire and Farnham and Surrey Heath Place Wound Management Formulary



Original Authors: Jane Harry TVN, Jane Turner TVN, Simon Smith MO technician.	V.1 Date signed off: 05/21
V.2 Amended by Simon Smith. Added Purilon to replace discontinued Aquaform. Added KCI negative pressure items.	V.2 Date signed off 06/22
V.3 Added on introductory paragraph	V3 Signed off 12/22
Date for next review: 05/23. Review by: Jane Harry TVN, Jane Turner TVN, Simon Smith MO technician. .	No Changes made
Date for next review 05/25	

All products on this formulary should be ordered through the “Online Non-Prescription Ordering Service” (ONPOS)

The ONPOS platform provides our healthcare professionals with an easy-to-use system to order wound care products. There is no requirement for a doctor to sign a prescription. The service can deliver products within 1-2 working days direct to your location.

If you need further information such as login details or how to use the service, please contact [simon.smith22@nhs.net](mailto:simon.smith22@nhs.net).

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

## Wound care formulary

**FOR ALL DRESSINGS/PRODUCTS: Store away from direct sunlight at ambient temperature and humidity. Do not use in cases where there is a known sensitivity to any of the ingredients.**

### Absorbent Cellulose

Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Xupad®</b> <b>Manufacturer:</b> Richardson Healthcare	<ul style="list-style-type: none"> <li>Moderate to heavily exuding wounds</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used as a primary dressing.</li> <li>Not to be used under compression.</li> </ul>	5 days	10cm x 12cm 10cm x 20cm 20cm x 20cm 20cm x 40cm	25 £3.75 25 £4.25 15 £4.20 12 £4.80	<b>DN, PN, NH, RH, TVN</b> <b>DN, PN, NH, RH, TVN</b> <b>DN, PN, NH, RH, TVN</b> <b>DN, PN, NH, RH, TVN</b>

### Alginate/Hydrofibres and hemostats

Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>First line</b> <b>Aquacel Extra and Aquacel Ribbon</b> <b>Manufacturer:</b> Convatec	<ul style="list-style-type: none"> <li>For the management of ulcers (leg, pressure (categories II–IV) and diabetic ulcers)</li> <li>Surgical wounds, donor sites, abrasions and lacerations.</li> <li>First-degree and second-degree burns.</li> <li>Traumatic wounds, painful wounds and wounds that are prone to bleeding.</li> <li>Use ribbon for dressing sinus, tracking or undermining wounds</li> </ul>	<ul style="list-style-type: none"> <li>Should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components</li> </ul>	Up to 7 days	5cmx5cm 10cmx10cm 15cmx15cm 1cmx45cm ribbon 2cmx45cm ribbon	10 £10.50 10 £25.10 5 £23.60 5 £9.55 5 £12.90	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b> <b>DN, PN, TVN</b> <b>DN, PN, NH, TVN</b> <b>DN, PN, TVN</b>
<b>Kaltostat®</b> <b>Manufacturer:</b> Convatec	<ul style="list-style-type: none"> <li>Acts as a haemostat</li> <li>Moderate to heavily exuding wounds and the local management of bleeding minor wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Should not be used on individuals who are sensitive or have had an allergic reaction to the dressing or its components.</li> <li>Not indicated for third-degree burns or use as a surgical sponge.</li> <li>Not indicated for heavily bleeding wounds.</li> </ul>	N/A	5cmx5cm 7.5cmx12cm	10 £10.00 10 £21.90	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

Antimicrobials *must be reviewed after two weeks*						
Silver Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Acticoat Flex 3®</b> <b>Manufacturer:</b> Smith and Nephew	<ul style="list-style-type: none"> <li>To be used under guidance from TVN, podiatry or vascular specialists.</li> <li>To reduce bioburden.</li> <li>Localised infection.</li> <li>Spreading infection.</li> <li>Systemic infection.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time.</li> <li>If wound is not showing signs of infection.</li> <li>Known sensitivity to silver.</li> <li>Patients undergoing MRI scan.</li> <li>Prior to radiotherapy.</li> <li>Where bioburden is not an issues</li> </ul>	3 days	5cm x 5cm 10cm x 10cm 10cmx20cm	5 £18.65 5 £45.55 12 £170.88	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b> <b>TVN only</b>
<b>Acticoat Flex 7®</b> <b>Manufacturer:</b> Smith and Nephew	<ul style="list-style-type: none"> <li>To be used under guidance from TVN, podiatry or vascular specialists.</li> <li>To reduce bioburden.</li> <li>Localised infection.</li> <li>Spreading infection.</li> <li>Systemic infection.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time, unless discussed with TVN</li> <li>If wound is not showing signs of infection.</li> <li>Known sensitivity to silver.</li> <li>Patients undergoing MRI scan.</li> <li>Prior to radiotherapy.</li> <li>Where bioburden is not an issues</li> </ul>	Up to 7 days	5cm x 5cm 10cm x 12.5cm 15cm x 15cm	5 £32.40 5 £96.55 5 £173.65	<b>TVN only</b> <b>TVN only</b> <b>TVN only</b>
<b>Aquacel Ag + Extra® and Ribbon</b> <b>Manufacturer:</b> ConvaTec	<ul style="list-style-type: none"> <li>To be used under guidance from TVN, not as first line antimicrobial.</li> <li>Moderately to highly exuding wounds.</li> <li>Ribbon to pack cavity wounds.</li> <li>To reduce bioburden.</li> <li>Localised infection.</li> <li>Spreading infection</li> <li>Systemic infection</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time.</li> <li>Dry wounds with minimal exudate.</li> <li>Known sensitivity to silver.</li> <li>Patients undergoing MRI scan.</li> <li>Prior to radiotherapy treatment.</li> <li>Where bioburden is not an issue.</li> </ul>	Up to 7 days	5cm x 5cm 10cm x 10cm 15cm x 15cm 1x45cm (ribbon) 2x45cm (ribbon)	10 £20.20 10 £48.20 5 £45.40 5 £16.20 5 £24.25	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b> <b>TVN only</b> <b>DN, PN, TVN</b> <b>DN, PN, TVN</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

DACC Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>First line</b> <b>Cutimed Sorbact Swabs®</b> <b>Manufacturer:</b> BSN	<ul style="list-style-type: none"> <li>Dressing for contaminated, colonised or infected wounds.</li> <li>For deep, superficial, traumatic, postoperative or dehisced wounds.</li> <li>For ulcers (venous, arterial, diabetic, pressure).</li> <li>Suitable for fungal infections in the groin, skin folds or between digits.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use in combination with ointments and creams as the binding effect is impaired</li> </ul>	N/A	4cm x 6cm 7cm x 9cm 2cm x 50cm (ribbon gauzes)	5 £8.70 5 £14.45 20 £84.80	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b> <b>TVN only</b>
<b>Cutimed Sorbact gel sheet®</b> <b>Manufacturer:</b> BSN	<ul style="list-style-type: none"> <li>Postoperative dehisced wounds, traumatic wounds and wounds occurring after excision of fistula.</li> <li>Chronic wounds such as venous, arterial, pressure or diabetic foot ulcers</li> </ul>	<ul style="list-style-type: none"> <li>Known hypersensitivity to the gel or its individual ingredients.</li> </ul>	N/A	7.5cm x 7.5cm 7.5cm x 15cm	10 £27.90 10 £47.10	<b>TVN only</b> <b>TVN only</b>
<b>Cutimed Sorbact Dressing pad®</b> <b>Manufacturer:</b> BSN	<ul style="list-style-type: none"> <li>Shallow, contaminated, colonised or infected wounds with low level of exudate.</li> <li>Also including superficial, traumatic, postoperative or dehisced wounds.</li> <li>Also for ulcers (venous, arterial, diabetic and pressure) and fungal infections</li> </ul>	<ul style="list-style-type: none"> <li>Do not use in combination with ointments and creams as the binding effect is impaired</li> </ul>	N/A	7cm x 9cm 10cm x 10cm 10cm x 20cm	5 £18.50 5 £28.95 20 £180.60	<b>TVN only</b> <b>TVN only</b> <b>TVN only</b>
Iodine Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Inadine®</b> <b>Manufacturer:</b> KCI/Acility	<ul style="list-style-type: none"> <li>To be used on superficial infected wounds e.g traumatic wounds, superficial burns.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time.</li> <li>Use with caution in people with thyroid disorders.</li> <li>Not to be used in moderate to wet wounds.</li> </ul>	3 days	5cm x 5cm 9.5cm x 9.5cm	25 £8.25 10 £4.90	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Iodoflex Paste®</b> <b>Manufacturer:</b> S&N	<ul style="list-style-type: none"> <li>For treatment of chronic exuding wounds.</li> <li>For infected wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time.</li> <li>To be used with caution in people who have thyroid disorders, pregnant or lactating women.</li> </ul>	7 days	5g 10g	5 £21.70 3 £25.98	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b>
<b>Honey Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Medihoney Apinate Dressing®</b> <b>Manufacturer:</b> Dermasciences	<ul style="list-style-type: none"> <li>Non adherent, non-absorbent, protease modulating matrix for use on infected wounds.</li> <li>Contains calcium and honey.</li> </ul>	<ul style="list-style-type: none"> <li>Honey should not be used on intact skin.</li> <li>Be aware that honey will increase exudate levels.</li> <li>Consider use of barrier skin protector to peri wound when using all Medihoney® dressings.</li> <li>Do not use honey dressings on wounds below the ankle without consulting a TVN/podiatrist.</li> </ul>	7 days	5cm x 5cm 10cm x 10cm 2cm x 30cm	10 £20.40 5 £17.30 5 £21.40	<b>DN, PN, NH, TVN</b> <b>DN, PN, TVN</b> <b>TVN only</b>
<b>Medihoney® Antibacterial Medical Honey Wound gel</b> <b>Manufacturer:</b> Dermasciences	<ul style="list-style-type: none"> <li>For use on infected wounds that require debridement of non-viable tissue.</li> </ul>	<ul style="list-style-type: none"> <li>Honey should not be used on intact skin.</li> <li>Be aware that honey will increase exudate levels.</li> <li>Consider use of barrier skin protector to peri wound when using all Medihoney® dressings.</li> <li>Do not use honey dressings on wounds below the ankle without consulting a TVN/podiatrist.</li> </ul>	7 days	20g	5 £20.15	<b>DN, PN, TVN</b>
<b>Medihoney® HCS Adhesive</b> <b>Manufacturer:</b> Dermasciences	<ul style="list-style-type: none"> <li>As above but with adhesive border.</li> </ul>	<ul style="list-style-type: none"> <li>Honey should not be used on intact skin.</li> <li>Be aware that honey will increase exudate levels</li> <li>Consider use of barrier skin protector to peri wound when using all Medihoney® dressings.</li> <li>Do not use honey dressings on wounds below the ankle without consulting a TVN/podiatrist.</li> </ul>	7 days	7.2cm x 7.2cm 11.5cm x 11.5cm	10 £31.20 10 £59.00	<b>TVN only</b> <b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

PHMB Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Suprasorb X + PHMB®</b> <b>Manufacturer:</b> L&R	<ul style="list-style-type: none"> <li>For use in light to moderately exuding wounds, requiring antimicrobials.</li> <li>For infected wounds hydrobalance dressing with antimicrobial.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used on full thickness burns, tendon, ligament or bone exposure.</li> </ul>	7 days	5cm x5cm 9cm x 9cm	5 £13.30 5 £26.50	<b>TVN only</b> <b>TVN only</b>
Octenilin Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Octenilin®</b> <b>Manufacturer:</b> Schulke	<ul style="list-style-type: none"> <li>Suitable for cleansing, moistening and decontamination of chronic skin wounds and to support the natural healing process.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on individuals who have known allergic reactions to any of the contents.</li> <li>Not to be used on hyaline cartilage.</li> </ul>	3-5 minutes irrigation	350ml 20ml wound gel	1 £4.60 1 £4.78	<b>DN, PN, TVN</b> <b>TVN only</b>
Bandages						
K-Four® Compression bandages.	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>K-Soft®</b> <b>Manufacturer:</b> Urgo	<ul style="list-style-type: none"> <li>Must be used under all leg bandaging</li> <li>Redistributes pressure</li> <li>1<sup>st</sup> layer of the k-four multilayer compression system</li> </ul>	<ul style="list-style-type: none"> <li>Allergy to lanolin</li> </ul>	7 days	10cmx3.5m 10cmx4.5m	1 £0.47 1 £0.59	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b>
<b>K-Lite®</b> <b>Manufacturer:</b> Urgo	<ul style="list-style-type: none"> <li>Light support bandage to aid absorbency</li> <li>2<sup>nd</sup> layer of the K-Four multilayer compression system</li> </ul>	<ul style="list-style-type: none"> <li>Known allergy to any of the components.</li> <li>Not suitable for application of pressure on its own, and so should be used in combination with other bandages in the treatment of venous leg ulcers</li> </ul>	7 days	10cmx4.5m 10cmx5.25m	1 £1.04 1 £1.19	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b>
<b>K-Plus®</b> <b>Manufacturer:</b> Urgo	<ul style="list-style-type: none"> <li>Provides 14-17mmHg at the ankle. Apply in figure of eight with 50% overlap and 50% stretch.</li> <li>3<sup>rd</sup> layer of the k-four multilayer compression system.</li> <li>Only for ankles under 25cm</li> <li>Only to be applied in a figure of 8 by a member of staff who have deemed competent.</li> </ul>	<ul style="list-style-type: none"> <li>Compression bandage systems are not recommended for patients with arterial disease</li> <li>Should not be used on patients with ABPI &lt;0.8 &gt;1.3</li> <li>Caution in diabetes, rheumatoid arthritis, congestive cardiac failure or peripheral neuropath</li> </ul>	7 days	10cmx8.7m 10cmx10.25m	1 £2.37 1 £2.74	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Ko-Flex®</b> <b>Manufacturer:</b> Urgo	<ul style="list-style-type: none"> <li>Provides 17-23mmHg at the ankle Apply in a spiral with 50% overlap &amp; 50% stretch</li> <li>4<sup>th</sup> layer of the k-four multilayer compression system</li> </ul>	<ul style="list-style-type: none"> <li>Patients with an allergy to latex</li> <li>No prior Doppler/Duplex assessment</li> <li>Compression bandage systems are not recommended for patients with arterial disease.</li> <li>Should not be used on patients with ABPI &lt;0.8 &gt;1.3</li> <li>Caution in diabetes, rheumatoid arthritis, congestive cardiac failure or peripheral neuropath</li> </ul>	7 days	10cmx6m 10cmx7m	1 £3.10 1 £3.60	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b>
<b>K-Three C®</b> <b>Manufacturer:</b> Urgo	<ul style="list-style-type: none"> <li>Type 3c high compression bandage. Provides 25–35mmHg and is a component of the K-Four 25–30cm, delivering the recommended pressure of 40mmHg when applied to larger ankles</li> <li>Only for ankles over 25cm</li> <li>Only to be applied in a spiral by a member of staff who have deemed competent.</li> </ul>	<ul style="list-style-type: none"> <li>Compression bandage systems are not recommended for patients with arterial disease.</li> <li>Should not be used on patients with ABPI &lt;0.8 &gt;1.3</li> <li>Caution in diabetes, rheumatoid arthritis, congestive cardiac failure or peripheral neuropath</li> </ul>	7 days	10cmx3m	1 £2.95	<b>DN, PN, TVN</b>
<b>Compression bandages.</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Urgo K-Two®</b> <b>Manufacturer:</b> Urgo	<ul style="list-style-type: none"> <li>Provides up to 40mmHg at ankle (full compression)</li> <li>Reduced compression kit provides 20mmHg at the ankle.</li> <li>Select size according to ankle circumference</li> <li>Latex free is available</li> </ul>	<ul style="list-style-type: none"> <li>Arterial disease</li> <li>Should not be used on patients with ABPI &lt;0.8 &gt;1.3</li> <li>Caution in diabetes, rheumatoid arthritis, congestive cardiac failure or peripheral neuropath</li> </ul>	7 days	18cm-25cm 25cm-32cm 18-25cm Latex Free 25-32cm Latex Free  <b>Reduced</b> 18cm-25cm 25cm-32cm 18-25cm Latex Free 25-32cm Latex Free	1 £7.88 1 £8.57 1 £9.13 1 £9.98  1 £8.87 1 £9.69 1 £9.13 1 £9.98	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b> <b>TVN</b> <b>TVN</b>  <b>DN, PN, TVN</b> <b>DN, PN, TVN</b> <b>TVN</b> <b>TVN</b>
<b>Actico®</b> <b>Manufacturer:</b> Lohmann & Rauscher	<ul style="list-style-type: none"> <li>Venous ulcers</li> <li>Lymphoedema</li> </ul>	<ul style="list-style-type: none"> <li>Caution if ankle circumference &lt;18cm with padding</li> <li>Should not be used on patients with ABPI &lt;0.8 &gt;1.3</li> <li>Caution in diabetes, rheumatoid arthritis, congestive cardiac failure or peripheral neuropath</li> </ul>	7 days	10cmx6m 8cmx6m 12cmx6m	1 £3.52 1 £3.39 1 £4.49	<b>DN, PN, TVN</b> <b>TVN only</b> <b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

Medicated bandages	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Viscopaste®</b> Manufacturer: S&N	<ul style="list-style-type: none"> <li>Chronic eczema or dermatitis</li> <li>Under graduated compression</li> </ul>	<ul style="list-style-type: none"> <li>Known sensitivity or allergy to any of the components.</li> </ul>	7 days	7.5cmx6m	1 £3.84	DN, PN, TVN
<b>Ichtopaste®</b> Manufacturer: Smith and Nephew	<ul style="list-style-type: none"> <li>Chronic eczema or dermatitis</li> <li>Under graduated compression</li> </ul>	<ul style="list-style-type: none"> <li>Known sensitivity or allergy to any of the components.</li> <li>May increase absorption of topical steroids, anaesthetics, retinoids</li> </ul>	7 days	7.5cmx6m	1 £3.87	TVN only
Retention bandages.	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>K-Band®</b> Manufacturer: Urgo	<ul style="list-style-type: none"> <li>Dressing retention</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	7 days	5cmx4m	1 £0.21	DN, PN, NH
Tubular bandages.	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Actifast®</b> Manufacturer: Lohmann & Rauscher	<ul style="list-style-type: none"> <li>Dressing retention</li> </ul>	<ul style="list-style-type: none"> <li>Not for routine use over bandaging.</li> <li>Not for routine use under compression.</li> </ul>	7 days	Beige 1m Blue 3m Blue 5m Yellow 3m Yellow 5m	1 £1.89 1 £2.20 1 £3.64 1 £3.61 1 £5.88	DN, PN, NH, TVN DN, PN, NH, RH, TVN DN, PN, NH, TVN DN, PN, NH, RH, TVN DN, PN, NH, TVN
<b>Comfinette</b> Manufacturer: Vernacare	<ul style="list-style-type: none"> <li>Can be used to hold dressings in place.</li> <li>Can be used for patients with known sensitivities to K Soft.</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	7 days	12-20m 78-20m	1 £4.16 1 £5.57	DN, PN, NH, TVN DN, PN, NH, TVN
Barrier cream/Skin Protectants						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Cavilon No Sting Barrier Film®</b> Manufacturer: 3M	<ul style="list-style-type: none"> <li>Primary barrier film</li> <li>Can be used on broken skin</li> <li>Barrier against bodily fluids.</li> <li>Protection around stoma sites and peri wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Allow product to dry completely before application of pads/clothes/dressings</li> </ul>	N/A	1ml 3ml	1 £4.09 1 £6.62	DN, PN, NH, TVN DN, PN, NH, TVN
<b>Cavilon Durable Barrier Cream®</b> Manufacturer: 3M	<ul style="list-style-type: none"> <li>To protect against bodily fluids and moisturised skin.</li> <li>Can be used on broken skin.</li> <li>For prevention of moisture associated skin damage</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on infected skin.</li> <li>Only to be used on the sacrum, buttocks, gluteal folds, perianal area and perineum.</li> </ul>	N/A	28g	20 £3.32	DN, PN, NH, TVN

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).



<b>Zerolon®</b> <b>Manufacturer:</b> Thorton & Ross	<ul style="list-style-type: none"> <li>To moisturise and protect intact, irritated, or inflamed skin.</li> <li>For use in the management of incontinence-associated dermatitis and Suitable for use with incontinence pads.</li> <li>To protect against bodily fluids such as urine, faeces, and exudate.</li> </ul>	<ul style="list-style-type: none"> <li>Hypersensitivity to any of the ingredients.</li> </ul>	N/A	28g	1 £2.33	<b>DN, PN, NH, TVN</b>
<b>Proshield Plus Protectant®</b> <b>Manufacturer:</b> Smith and Nephew	<ul style="list-style-type: none"> <li>Intact and broken skin and skin damaged from incontinence</li> <li>Can also be used to moisturise several dry skin.</li> </ul>	<ul style="list-style-type: none"> <li>Only to be used on the sacrum, buttocks, gluteal folds, perianal area and perineum.</li> </ul>	N/A	115g	1 £9.94	<b>DN, PN, NH, TVN</b>
<b>Cavilon Advance®</b> <b>Manufacturer:</b> 3M	<ul style="list-style-type: none"> <li>Protection of intact or damaged skin from bodily fluids.</li> <li>Prevention and reversal of all forms of moisture associated skin damage.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on full-thickness wounds, around the eyes or on infected skin.</li> <li>Extremely flammable until completely dried on the skin.</li> </ul>	72 hours	2.7ml	1 £10.59	<b>TVN only</b>
<b>Lifteez medical adhesive remover</b> <b>Manufacturer:</b> Medicare plus	<ul style="list-style-type: none"> <li>For ease of removal of adhesive dressings, pouches and medical devices and to minimise the potential for skin stripping and pain on removal.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on infected areas of skin or if any irritation occurs.</li> <li>Avoid use on delicate or sensitive areas of skin.</li> <li>Only use as directed.</li> </ul>	N/A	50mls	1 £5.60	<b>TVN only</b>
<b>Medihoney Barrier Cream®</b> <b>Manufacuturer:</b> Medicare plus	<ul style="list-style-type: none"> <li>Protects at-risk skin from breakdown.</li> <li>Protects wound margins from exudate damage.</li> <li>Protects skin from damage caused by friction or shear.</li> <li>Treatment of incontinence-associated dermatitis and atopic eczema.</li> <li>Treats blistered skin.</li> <li>Prevention and treatment of maceration and excoriation.</li> <li>Helps maintain the skin pH.</li> <li>Suitable for all skin types.</li> </ul>	<ul style="list-style-type: none"> <li>Patch test before first use.</li> </ul>	N/A	50g	1 £4.89	<b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

Charcoal Dressings						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Clinisorb®</b> <b>Manufacturer:</b> Clinimed	<ul style="list-style-type: none"> <li>Management of malodorous wounds such as fungating wounds, pressure ulcers, leg ulcers and diabetic foot ulcers.</li> </ul>	<ul style="list-style-type: none"> <li>Store away from direct sunlight at ambient temperature and humidity.</li> <li>Do not use in cases where there is a known sensitivity to any ingredients.</li> </ul>	N/A	10cmx10cm 10cmx20cm	10 £19.90 10 £26.40	<b>TVN only</b> <b>TVN only</b>
Clear Acrylic						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Tegaderm Absorbent Clear Acrylic®</b> <b>Manufacturer:</b> 3M	<ul style="list-style-type: none"> <li>Skin tears</li> <li>Low to moderately exuding acute wounds</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	N/A	7.6cmx9.5cm 11.1cmx12.7cm	5 £16.15 5 £20.90	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b>
Contact layers						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Adaptic Touch®</b> <b>Manufacturer:</b> KCL	<ul style="list-style-type: none"> <li>For ulcer (Leg, pressure and diabetic ulcers)</li> <li>For donor sites</li> <li>1 st &amp; 2nd degree burns</li> <li>Traumatic (skin tears) and surgical wounds</li> <li>Allows free passage of exudate</li> <li>Requires secondary dressing</li> </ul>	<ul style="list-style-type: none"> <li>Known sensitivity to any of the dressing components.</li> </ul>	7 days	5cmx7.6cm 7.6cmx11cm 12.7cmx15cm 20cmx30cm	10 £11.30 10 £22.50 10 £46.50 5 £62.50	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b> <b>TVN only</b> <b>TVN only</b>
<b>N-A Ultra®</b> <b>Manufacturer:</b> KCL	<ul style="list-style-type: none"> <li>Ulcers (leg and pressure)</li> <li>Burns</li> <li>Other granulating wounds</li> <li>Allows free passage of exudate</li> <li>Requires secondary dressing</li> </ul>	<ul style="list-style-type: none"> <li>Do not use if allergic to silicone.</li> </ul>	Up to 7 days	9.5cmx9.5cm 9.5cmx19cm	40 £13.60 25 £16.00	<b>DN, PN, NH, RH, TVN</b> <b>DN, PN, NH, RH, TVN</b>
<b>Urgostart Plus Pad®</b> <b>Manufacturer:</b> Urgo medical	<ul style="list-style-type: none"> <li>For all stages (from desloughing to complete healing) of exuding wounds including leg ulcers, diabetic foot ulcers, pressure ulcers, and long-standing acute wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on heavily bleeding wounds, cancerous wounds, wounds that may reveal a deep abscess, infected wounds.</li> <li>Where there is a known sensitivity to the dressing.</li> </ul>	Up to 7 days	10cmx10cm	10 £50.70	<b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Urgostart® contact Manufacturer:</b> Urgo medical	<ul style="list-style-type: none"> <li>Requires secondary dressing</li> <li>Chronic wounds</li> <li>Cavity wounds</li> <li>Can be used under compression</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on infected or critically colonised wounds</li> <li>Cancerous wounds</li> <li>Fistulas, which may reveal a deep abscess</li> <li>Known sensitivity to the dressing.</li> </ul>	Up to 7 days	10cmx10cm	10 £43.30	<b>TVN only</b>
<b>Urgotul® Manufacturer:</b> Urgo medical	<ul style="list-style-type: none"> <li>Leg and pressure ulcers</li> <li>1st &amp; 2nd degree burns</li> <li>Traumatic (skin tears), donor site and surgical wounds</li> <li>Allows free passage of exudate</li> <li>Requires secondary dressing</li> </ul>	<ul style="list-style-type: none"> <li>Known sensitivity to any components.</li> <li>Product contains hydrocolloid and petroleum jelly particles</li> </ul>	Up to 7 days	5cm5cm 10cmx10cm 15cmx20cm 20cmx30cm	10 £15.70 10 £31.30 10 £88.60 5 £71.25	<b>All TVN only</b>
<b>Dressing packs</b>						
<b>Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Nurse it dressing pack</b>			N/A	Small/Medium Medium/Large	10 £5.55 10 £5.55	<b>DN, PN, NH, TVN DN, PN, NH, TVN</b>
<b>Emollients</b>						
<b>Items</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Zerocream® Cream Manufacturer:</b> Zeroderma	<ul style="list-style-type: none"> <li>Eczema</li> <li>Psoriasis</li> <li>Dry skin conditions</li> </ul>	<ul style="list-style-type: none"> <li>Known sensitivity or allergy to any components.</li> <li>Potential fire hazard with paraffin based emollients in contact with dressings or clothing as they maybe ignited by naked flame</li> </ul>	N/A	500g	1 £4.08	<b>DN, PN, NH, TVN</b>
<b>Zeroderm® ointment Manufacturer:</b> Zeroderma	<ul style="list-style-type: none"> <li>Eczema</li> <li>Psoriasis</li> <li>Dry skin conditions</li> </ul>	<ul style="list-style-type: none"> <li>Known sensitivity or allergy to any components.</li> <li>Potential fire hazard with paraffin based emollients in contact with dressings or clothing as they maybe ignited by naked flame</li> </ul>	N/A	500g	1 £4.10	<b>DN, PN, NH, TVN</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

Films						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>ClearFilm®</b> <b>Manufacturer:</b> Richardson Healthcare	<ul style="list-style-type: none"> <li>Minor burns</li> <li>Protective cover and fixation of catheter sites</li> <li>Skin graft donor sites</li> <li>Clean closed surgical incisions</li> <li>Abrasions</li> <li>Blisters</li> <li>Secondary dressing</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	7 days	6cmx7cm 10cmx12cm	50 £12.50 10 £5.50	DN, PN, NH, TVN DN, PN, NH, TVN
Foams						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>First line</b> <b>Biatain</b> <b>Silicone®</b> <b>Manufacturer:</b> Coloplast	<ul style="list-style-type: none"> <li>Primary dressing does not require contact layer.</li> <li>Moderately exuding wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Don't use with oxidising solutions such as hypochlorite or hydrogen peroxide.</li> </ul>	Up to 7 days	7.5cm x 7.5cm 10cm x 10cm 12.5cm x 12.5cm 15cm x 15cm 17.5cm x 17.5cm	10 £15.10 10 £22.30 10 £27.20 10 £40.40 10 £53.70	DN, PN, NH, RH, TVN DN, PN, NH, RH, TVN DN, PN, NH, RH, TVN DN, PN, TVN DN, PN, TVN
<b>Aquacel Foam</b> <b>Adhesive®</b> <b>Manufacturer:</b> ConvaTec	<ul style="list-style-type: none"> <li>Primary dressing and does not require contact layer.</li> <li>Moderately exuding wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Dry wounds with minimal exudate.</li> <li>Do not moisten prior to application.</li> </ul>	Up to 7 days.	8cm x 8cm 8cm x 13cm 10cm x 10cm 10cm x 20cm 10cm x 30cm 12.5cm x 12.5cm 17.5cm x 17.5cm 21cm x 21cm 25cm x 30cm	10 £14.70 10 £23.50 10 £22.50 5 £13.70 5 £26.05 10 £27.90 10 £55.80 5 £40.85 5 £52.90	DN, PN, NH, TVN DN, PN, NH, TVN DN, PN, NH, TVN DN, PN, TVN DN, PN, TVN DN, PN, NH, TVN TVN only TVN only TVN only
<b>Aquacel Foam</b> <b>Non-Adhesive®</b> <b>Manufacturer:</b> ConvaTec	<ul style="list-style-type: none"> <li>Primary dressing and does not require contact layer.</li> <li>Moderately exuding wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Dry wounds with minimal exudate.</li> <li>Do not moisten prior to application.</li> </ul>	Up to 7 days.	5cm x 5cm 10cm x 10cm 15cm x 15cm 20cm x 20cm 15cm x 20cm	10 £14.10 10 £26.60 5 £22.40 5 £36.50 5 £30.60	DN, PN, NH, TVN DN, PN, NH, TVN DN, PN, TVN DN, PN, TVN DN, PN, TVN
<b>Biatain soft</b> <b>hold®</b> <b>Manufacturer:</b> Coloplast	<ul style="list-style-type: none"> <li>Primary dressing and does not require contact layer.</li> <li>Moderately exuding wounds.</li> <li>Can be used under compression therapy</li> </ul>	<ul style="list-style-type: none"> <li>Do no use with oxidising solutions such as hypochlorite or hydrogen peroxide.</li> </ul>	Up to 7 days	5cm x 7cm 10cm x 10cm 10cm x 20cm	5 £6.90 5 £13.60 5 £20.65	DN, PN, TVN DN, PN, TVN DN, PN, TVN

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Biatain Adhesive Sacral®</b> <b>Manufacturer:</b> Coloplast	<ul style="list-style-type: none"> <li>Primary dressing and does not require contact layer.</li> <li>Moderately exuding wounds.</li> <li>Can be useful when having issues with adhesion of other dressing to the sacral area,</li> </ul>	<ul style="list-style-type: none"> <li>Dry wounds with minimal exudate.</li> <li>Do no use with oxidising solutions such as hypochlorite or hydrogen peroxide.</li> </ul>	Up to 7 days	17cm x 17cm 23cm x 23cm	5 £12.20 5 £23.20	<b>TVN only</b> <b>TVN only</b>
<b>Biatain Non-Adhesive®</b> <b>Manufacturer:</b> Coloplast	<ul style="list-style-type: none"> <li>Primary dressing and does not require contact layer.</li> <li>Moderately exuding wounds.</li> <li>Can be used under compression therapy</li> </ul>	<ul style="list-style-type: none"> <li>Dry wounds with minimal exudate.</li> <li>Do no use with oxidising solutions such as hypochlorite or hydrogen peroxide</li> </ul>	Up to 7 days	5cm x 7cm 10cm x 10cm 15cm x 15cm 10cm x 20cm	10 £13.80 10 £25.00 5 £23.05 5 £20.65	<b>TVN only</b> <b>TVN only</b> <b>TVN only</b> <b>TVN only</b>
Hydrocolloids						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Comfeel plus transparent®</b> <b>Manufacturer:</b> Coloplast	<ul style="list-style-type: none"> <li>Partial and full thickness wounds</li> <li>Wounds with no, or low exudate</li> <li>Draws splinters</li> <li>Supports autolytic debridement (to be used for this with TVN input/discussion)</li> <li>Maintains a moist environment</li> </ul>	<ul style="list-style-type: none"> <li>Highly exuding wounds</li> <li>Presence of infection</li> <li>Diabetic foot ulcers</li> </ul>	Up to 7 days	5cmx7cm 10cmx10cm 15cmx15cm	10 £6.90 10 £13.20 5 £17.20	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b> <b>TVN only</b>
<b>Duoderm Extra Thin®</b> <b>Manufacturer:</b> ConvTec	<ul style="list-style-type: none"> <li>Partial and full thickness wounds</li> <li>Wounds with no or low exudate</li> <li>Draws splinters</li> <li>Supports autolytic debridement (to be used for this with TVN input/discussion)</li> <li>Maintains a moist environment</li> </ul>	<ul style="list-style-type: none"> <li>Highly exuding wounds</li> <li>Presence of infection</li> <li>Diabetic foot ulcers</li> <li>Contains gelatine derived from pork. Consider patients with religious or ethical objections</li> </ul>	Up to 7 days	10cmx10cm 15cmx15cm	10 £13.90 10 £30.00	<b>DN, TVN</b> <b>DN, TVN</b>
Hydrogels						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Purilon®</b> <b>Manufacturer:</b> Coloplast	<ul style="list-style-type: none"> <li>Primarily intended for treatment of necrotic and sloughy wounds, e.g. leg ulcers, pressure ulcers and non-infected diabetic foot ulcers.</li> <li>First and second-degree burns.</li> <li>Can be used throughout the healing process to provide a moist healing environment.</li> </ul>	<ul style="list-style-type: none"> <li>Not recommended on wet wounds</li> </ul>	N/A	8g	1 £1.85	<b>DN, PN, TVN</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Kerralite Cool Non-Adhesive®</b> <b>Manufacturer:</b> Crawford	<ul style="list-style-type: none"> <li>Chronic wounds and painful wounds,</li> <li>Skin conditions such as leg ulcers, radiation therapy damage, skin tears, burns and scalds.</li> <li>May be used on low-exuding and non-exuding wounds to assist in autolytic debridement by hydration of necrotic and sloughy tissue and for absorption of exudate</li> <li>This product needs to be cut to the size of the wound.</li> </ul>	<ul style="list-style-type: none"> <li>Full-thickness wounds</li> <li>Heavily bleeding wounds</li> <li>Third-degree burns</li> <li>A covering for deep, narrow cavities or sinuses</li> <li>Not to be used on the lower leg, or foot without consultation with TVN/Podiatry</li> </ul>	N/A	6cmx6cm 12cmx8.5cm 18cmx12.5cm	5 £8.95 5 £13.20 5 £19.10	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b> <b>TVN only</b>
Miscellaneous						
Cleaning Solutions	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>SAL-E-Pods</b> <b>Manufacturer:</b> Ennogen	<ul style="list-style-type: none"> <li>For topical irrigation and cleansing of wounds.</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	N/A	20ml	1 £4.80	<b>DN, PN, NH, RH, TVN</b>
Debridement	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Debrisoft®</b> <b>Manufacturer:</b> L&R	<ul style="list-style-type: none"> <li>For removal of superficial slough and debris, including biofilm, from wounds and skin during wound bed preparation</li> <li>Can also be used at dressing change to aid assessment.</li> </ul>	<ul style="list-style-type: none"> <li>Prior treatment with autolytic debridement is advised if very stubborn slough or hard necrosis is present.</li> <li>Not to be used in conjunction with emollients as will not be as effective.</li> </ul>	N/A	10cmx10cm	5 £33.45	<b>DN, PN, TVN</b>
Gauze and swabs	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Non-woven swabs</b> <b>Manufacturer:</b> Medicare Plus		<ul style="list-style-type: none"> <li>Not a leave on product</li> </ul>	N/A	10cmx10cm	100 £0.85	<b>DN, PN, NH, TVN</b>
<b>Non-woven swabs sterile</b> <b>Manufacturer:</b> Medicare Plus		<ul style="list-style-type: none"> <li>Not a leave on product</li> </ul>	N/A	7.5cmx7.5cm	5 £5.40	<b>DN, PN, NH, RH, TVN</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

Items	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Acti-Glide®</b> <b>Manufacturer:</b> L&R	<ul style="list-style-type: none"> <li>Aids application of compression hosiery (below-knee and thigh length) and arm sleeves.</li> </ul>	<ul style="list-style-type: none"> <li>Known sensitivity to the fabric.</li> </ul>	N/A	One size	1 £15.22	<b>TVN only</b>
<b>Leukostrip®</b> <b>Manufacturer:</b> Smith and Newpew	<ul style="list-style-type: none"> <li>For closure of minor wounds and minor operations in a primary care setting.</li> </ul>	<ul style="list-style-type: none"> <li>Skin tears</li> </ul>	3 days	6.4cmx76mm	10 packs of 3 strips £6.64	<b>DN, PN, TVN</b>
<b>LimBo Waterproof Protector®</b> <b>Manufacturer:</b> Thesis	<ul style="list-style-type: none"> <li>Used to keep leg or foot dressings, PICC lines, midlines, and plaster casts dry during bathing or showering.</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	N/A	Small build Slim build Normal-half leg Normal-short leg	1 £10.20 1 £10.20 1 £10.20 1 £10.20	<b>All TVN only</b>
<b>Simple dressings</b>						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>Clearpore®</b> <b>Manufacturer:</b> Richardson	<ul style="list-style-type: none"> <li>Low to moderately exuding wounds.</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	5 days	6cmx7cm 10cmx10cm 10cmx15cm	60 £7.20 50 £10.00 50 £12.00	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b>
<b>Softpore®</b> <b>Manufacturer:</b> Richardson	<ul style="list-style-type: none"> <li>Low to moderately exuding wounds.</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	5 days	6cmx7cm 10cmx10cm 10cmx20cm	60 £3.60 50 £6.50 30 £10.50	<b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b> <b>DN, PN, NH, TVN</b>
<b>Superabsorbants</b>						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>First line</b> <b>Vliwasorb Pro®</b> <b>Manufacturer:</b> L&R	<ul style="list-style-type: none"> <li>Heavily exuding wounds such as leg ulcers, pressure ulcers, diabetic foot ulcers, oncology wounds, partial-thickness burns, postoperative wounds, traumatic wounds and skin graft donor sites.</li> <li>Can be used as a primary dressing</li> <li>Also suitable for lymphorrhoea.</li> <li>Can be used under compression.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on lightly or non-exuding wounds.</li> <li>Not for use in cavity wounds as the dressing will expand on absorption of exudate.</li> <li>Do not cut.</li> </ul>	N/A	12.5cm x 12.5cm 12.5cm x 22.5cm 22cm x 22cm	10 £9.20 10 £10.80 10 £19.40	<b>DN, PN, TVN</b> <b>DN, PN, TVN</b> <b>DN, PN, TVN</b>
<b>Biatain Super Adhesive®</b> <b>Manufacturer:</b> Coloplast	<ul style="list-style-type: none"> <li>Highly exuding wounds</li> <li>Primary dressing.</li> <li>Infected wounds.</li> <li>Can be used under compression</li> </ul>	<ul style="list-style-type: none"> <li>Dry wounds</li> <li>Do not cut</li> </ul>	5-7 days	12.5cm x 12.5cm 15cm x 15cm 20cm x 20cm	10 £37.30 10 £45.00 10 £70.20	<b>TVN only</b> <b>TVN only</b> <b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Cutimed Sorbion Sachet S®</b> <b>Manufacturer:</b> BSN	<ul style="list-style-type: none"> <li>For wound cleansing and slough reduction of wounds with moderate or high exudate levels.</li> <li>Can be used as primary dressing</li> </ul>	<ul style="list-style-type: none"> <li>Avoid using on wounds with very little or no fluids (risk of drying and adhesion).</li> <li>Only apply to tunnel cavities under appropriate monitoring conditions as the product expands as it absorbs fluid.</li> </ul>	N/A	7.5cm x 7.5cm 10cm x 10cm 20 cm x 10cm 20cm x 20cm 30cm x 20cm	10 £9.90 10 £15.50 10 £20.70 10 £36.20 10 £51.80	<b>TVN only</b> <b>TVN only</b> <b>TVN only</b> <b>TVN only</b> <b>TVN only</b>
<b>Cutimed Sorbion Sachet Multi Star®</b> <b>Manufacturer:</b> BSN	<ul style="list-style-type: none"> <li>Management of moderate to highly exuding wounds</li> <li>Useful for hard to dressed areas such as axillae</li> <li>Can be used as primary dressing</li> </ul>	<ul style="list-style-type: none"> <li>Avoid using on wounds with very little or no fluids (risk of drying and adhesion).</li> <li>Only apply to tunnel cavities under appropriate monitoring conditions as the product expands as it absorbs fluid.</li> </ul>	N/A	8cm x 8cm 14cm x 14cm	10 £31.10 10 £50.80	<b>TVN only</b> <b>TVN only</b>
<b>Cutimed Sorbion Sachet XL®</b> <b>Manufacturer:</b> BSN	<ul style="list-style-type: none"> <li>Management of moderate to highly exuding wounds</li> </ul>	<ul style="list-style-type: none"> <li>Avoid using on wounds with very little or no fluids (risk of drying and adhesion).</li> <li>Only apply to tunnel cavities under appropriate monitoring conditions as the product expands as it absorbs fluid.</li> </ul>	N/A	45cm x 25cm	5 £77.80	<b>TVN only</b>
<b>Kerramax care®</b> <b>Manufacturer:</b> Crawford	<ul style="list-style-type: none"> <li>Highly exuding wounds</li> <li>Primary dressing</li> </ul>	<ul style="list-style-type: none"> <li>Dry wounds</li> <li>Avoid using on wounds with very little or no fluids (risk of drying and adhesion).</li> <li>Only apply to tunnel cavities under appropriate monitoring conditions as the product expands as it absorbs fluid</li> </ul>	N/A	10cm x 10cm 10cm x 22cm 20cm x 22cm 20cm x 30cm	10 £13.10 10 £17.10 10 £30.20 10 £34.50	<b>TVN only</b> <b>TVN only</b> <b>TVN only</b> <b>TVN only</b>
<b>Surgical Tape</b>						
<b>Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>First line</b> <b>Clinipore®</b> <b>Manufacturer:</b> Clinisupplies	<ul style="list-style-type: none"> <li>Tape to secure bandaging.</li> </ul>	<ul style="list-style-type: none"> <li>Allergy to tape adhesive</li> </ul>	N/A	2.5cmx5m	1 £0.60	<b>DN, PN, NH, RH, TVN</b>
<b>Second line</b> <b>Chemifix®</b> <b>Manufacturer:</b> Medicareplus	<ul style="list-style-type: none"> <li>Use with those patients that have allergy to Clinipore®, or those patients that require tape in direct contact with skin over prolonged periods of time.</li> </ul>	<ul style="list-style-type: none"> <li>Allergy to tape adhesive</li> </ul>	N/A	5cmx5m	1 £1.25	<b>DN, PN, TVN</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).



<b>Total Negative Pressure</b> <b>N.B - Only to be used in consultation with TVN services, for assessment of suitability and then guidance on product requirements and usage. Please contact vcl.tissueviability@nhs.net</b>						
Dressings	Indications	Contraindications	Maximum wear time	Size	Cost per unit/pack size	Available to
<b>ActiV.A.C. Canister®</b> Manufacturer: 3M+KCI	<ul style="list-style-type: none"> <li>For use with the ActiV.A.C. Therapy Unit.</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	N/A	300ml	1 £29.71	<b>TVN only</b>
<b>Kerlix AMD® Gauze</b> Manufacturer: H&R Healthcare	<ul style="list-style-type: none"> <li>For use as a primary wound filler under negative pressure.</li> </ul>	<ul style="list-style-type: none"> <li>Can be used in conjunction with prescribed therapies for the treatment of infection.</li> <li>Not intended as a primary treatment for infection or third-degree burns.</li> <li>Do not use on patients with known sensitivity to polyhexamethylene biguanide (PHMB).</li> </ul>	N/A	11.4cmx3.7m roll	1 £1.64	<b>TVN only</b>
<b>Pico 7®</b> Manufacturer: Smith and Nephew	<ul style="list-style-type: none"> <li>Chronic, acute, traumatic, subacute and dehiscent wounds</li> <li>Partial-thickness burns</li> <li>Ulcers (such as diabetic or pressure)</li> <li>Flaps and grafts, and surgically closed incision sites.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on patients with malignancy in the wound bed, or wound margins (except in palliative care to improve quality of life)</li> <li>Previously confirmed or untreated osteomyelitis</li> <li>Non-enteric and unexplored fistulas</li> <li>Necrotic tissue with eschar present</li> <li>Exposed arteries, veins, nerves or organs</li> <li>Exposed anastomotic sites</li> <li>emergency airway aspiration</li> <li>Pleural, mediastinal or chest tube drainage; and surgical suction.</li> </ul>	N/A	10cmx20cm 15cmx15cm 15cmx20cm	1 £133.06 1 £132.29 1 £132.29	<b>TVN only</b> <b>TVN only</b> <b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Renasys G® with soft port dressings</b> <b>Manufacturer:</b> Smith and Nephew	<ul style="list-style-type: none"> <li>• A range of wound types including postoperative and dehisced wounds</li> <li>• Ulcers including pressure ulcers, diabetic/neuropathic ulcers and insufficiency ulcers</li> <li>• Traumatic wounds,</li> <li>• Skin flaps and grafts</li> <li>• Venous and explored fistulae.</li> </ul>	<ul style="list-style-type: none"> <li>• Do not use on patients with malignancy in the wound bed or wound margins (except in palliative care to improve quality of life)</li> <li>• Previously confirmed or untreated osteomyelitis</li> <li>• Non-enteric and unexplored fistulas</li> <li>• Necrotic tissue with eschar present</li> <li>• Exposed arteries, veins, nerves or organs</li> <li>• Exposed anastomotic sites</li> <li>• Emergency airway aspiration</li> <li>• Pleural, mediastinal or chest tube drainage; and surgical suction.</li> </ul>	N/A	Small Medium Large	1 £20.41 1 £25.60 1 £32.48	<b>TVN only</b> <b>TVN only</b> <b>TVN only</b>
<b>Renasys F with soft port dressing kit®</b> <b>Manufacturer:</b> Smith and Nephew	<ul style="list-style-type: none"> <li>• A range of wound types including postoperative and dehisced wounds</li> <li>• Ulcers including pressure ulcers, diabetic/neuropathic ulcers and insufficiency ulcers</li> <li>• Traumatic wounds,</li> <li>• Skin flaps and grafts</li> <li>• Venous and explored fistulae.</li> </ul>	<ul style="list-style-type: none"> <li>• Do not use on patients with malignancy in the wound bed or wound margins (except in palliative care to improve quality of life)</li> <li>• Previously confirmed or untreated osteomyelitis</li> <li>• Non-enteric and unexplored fistulas</li> <li>• Necrotic tissue with eschar present; exposed arteries, veins, nerves or organs.</li> <li>• Exposed anastomotic sites</li> <li>• Emergency airway aspiration.</li> <li>• Pleural, mediastinal or chest tube drainage; and surgical suction.</li> </ul>	N/A	Small Medium Large	1 £23.95 1 £27.83 1 £33.01	<b>All TVN only</b>
<b>Renasys Touch Canister®</b> <b>Manufacturer:</b> Smith and Nephew	<ul style="list-style-type: none"> <li>• Canister for fluid collection for Renasys Touch machine</li> </ul>	<ul style="list-style-type: none"> <li>• Do not use on patients with malignancy in the wound bed or wound margins (except in palliative care to improve quality of life)</li> <li>• Previously confirmed or untreated osteomyelitis</li> <li>• Non-enteric and unexplored fistulas</li> <li>• Necrotic tissue with eschar present</li> <li>• Exposed arteries, veins, nerves, or organs</li> <li>• Exposed anastomotic sites</li> <li>• emergency airway aspiration</li> </ul>	N/A	300ml	1 £29.04	<b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

		<ul style="list-style-type: none"> <li>Pleural, mediastinal or chest tube drainage; and surgical suction.</li> </ul>				
<b>Renasys Adhesive Gel Patch®</b> <b>Manufacturer:</b> Smith and Nephew	<ul style="list-style-type: none"> <li>For use with topical negative pressure</li> <li>To form an adhesive seal in support of existing dressings</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on patients with malignancy in the wound bed or wound margins (except in palliative care to improve quality of life)</li> <li>Previously confirmed or untreated osteomyelitis</li> <li>Non-enteric and unexplored fistulas</li> <li>Necrotic tissue with eschar present</li> <li>Exposed arteries, veins, nerves or organs</li> <li>Exposed anastomotic sites</li> <li>Emergency airway aspiration</li> <li>Pleural, mediastinal or chest tube drainage; and surgical suction.</li> </ul>	N/A	10x10cm	1 £3.56	<b>TVN only</b>
<b>T.R.A.C. Y Connector®</b> <b>Manufacturer:</b> 3M+KCI	<ul style="list-style-type: none"> <li>Compatible with all V.A.C. therapy dressings</li> </ul>	<ul style="list-style-type: none"> <li>SensaT.R.A.C. technology senses only one wound site.</li> <li>Avoid using a Y-Connector to connect wounds that are being optimally treated with different pressure settings.</li> <li>Not recommended for use with flaps or grafts.</li> </ul>	N/A	One size	1 £3.27	<b>TVN only</b>
<b>V.A.C. Drape®</b> <b>Manufacturer:</b> 3M+KCI	<ul style="list-style-type: none"> <li>Compatible with all V.A.C. therapy dressings.</li> </ul>	<ul style="list-style-type: none"> <li>Contains an acrylic adhesive coating that may present a risk of an adverse reaction in patients who are allergic or hypersensitive to acrylic adhesives.</li> </ul>		30.5x26cm	1 £9.70	<b>TVN only</b>
<b>V.A.C. Gel Strip®</b> <b>Manufacturer:</b> 3M+KCI	<ul style="list-style-type: none"> <li>For use with V.A.C. therapy to help seal negative pressure leaks.</li> </ul>	<ul style="list-style-type: none"> <li>None listed</li> </ul>	N/A	14x3cm	1 £3.93	<b>TVN only</b>
<b>V.A.C. Granufoam Bridge Dressing Kit®</b> <b>Manufacturer:</b> 3M+KCI	<ul style="list-style-type: none"> <li>V.A.C. GranuFoam can be used on most wound types: chronic, acute traumatic, sub-acute, partial-thickness burns, pressure ulcers, venous and diabetic foot ulcers, flaps and grafts.</li> <li>For use when it is necessary to place the SensaT.R.A.C.</li> </ul>	<ul style="list-style-type: none"> <li>For use only when it is necessary to place the SensaT.R.A.C. Pad away from over the wound site.</li> </ul>		37.5cm x 32.5cm x 25.5cm	5 £167.50	<b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

	<ul style="list-style-type: none"> <li>Pad away from over the wound site.</li> </ul>					
<b>V.A.C. Granufoam Dressing kits with SENSAT.R.A.C. Manufacturer: 3M+KCI</b>	<ul style="list-style-type: none"> <li>Can be used on most wound types: chronic, acute traumatic, sub-acute, partial-thickness burns, venous and diabetic foot ulcers, diabetic ulcers, flaps and grafts.</li> </ul>	<ul style="list-style-type: none"> <li>Do not place over exposed blood vessels or organs.</li> </ul>		Small Medium Large	1 £24.00 5 £120.00 1 £28.57 5 £139.15 1 £33.14 5 £165.70	<b>TVN only</b>
<b>V.A.C. Granufoam Silver Dressing kits with SENSAT.R.A.C. Manufacturer: 3M+KCI</b>	<ul style="list-style-type: none"> <li>Consider for use on infected wounds or wounds at risk of infection.</li> </ul>	<ul style="list-style-type: none"> <li>Sensitivity to silver; do not place V.A.C. dressing over exposed blood vessels or organs.</li> </ul>		10cmx7.5cm 18cmx12.5cm	5 £167.05 10 £380.40	<b>TVN only</b>
<b>V.A.C. Simplace EX Dressing Kit Manufacturer: 3M+KCI</b>	<ul style="list-style-type: none"> <li>Can be used on most wound types: chronic, acute traumatic, sub-acute, partial-thickness burns, pressure ulcers, venous leg ulcers, diabetic foot ulcers, flaps and grafts.</li> </ul>	<ul style="list-style-type: none"> <li>Exposed blood vessels or organs.</li> </ul>		Small Medium	5 £139.05 5 £159.90	<b>TVN only</b>
<b>V.A.C Whitefoam Dressing® Manufacturer: 3M+KCI</b>	<ul style="list-style-type: none"> <li>For wounds that contain areas of undermining, sinuses or tunnels. Also suitable for split-thickness or full-thickness skin grafts.</li> <li>Recommended for use in wounds where the growth of granulation tissue into the foam needs to be controlled.</li> </ul>	<ul style="list-style-type: none"> <li>Do not place V.A.C. dressing over exposed blood vessels or organs.</li> </ul>		Small	10 £111.30	<b>TVN only</b>
<b>V.A.C Whitefoam Dressing Kit® Manufacturer: 3M+KCI</b>	<ul style="list-style-type: none"> <li>Wounds that contain areas of undermining, sinuses or tunnels.</li> <li>Also suitable for split-thickness or full-thickness skin grafts.</li> <li>Recommended for use in wounds where the growth of granulation</li> </ul>	<ul style="list-style-type: none"> <li>Do not place V.A.C. dressing over exposed blood vessels or organs.</li> </ul>	N/A	Small	1 £27.09	<b>TVN only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

	tissue into the foam needs to be controlled.					
<b>Podiatry formulary</b>						
<b>Alginate</b>						
<b>Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Kaltostat®</b> <b>Manufacturer:</b> ConvaTec	<ul style="list-style-type: none"> <li>Acts as a haemostat</li> <li>Moderate to heavily exuding wounds and the local management of bleeding minor wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Should not be used on individuals who are sensitive or have had an allergic reaction to the dressing or its components.</li> <li>Not indicated for third-degree burns or use as a surgical sponge.</li> <li>Not indicated for heavily bleeding wounds.</li> </ul>	7 days	5cm x5cm	10 £10.00	<b>DN, PN, NH, TVN/Pod</b>
<b>Antimicrobials *must be reviewed after six weeks use on the foot*</b>						
<b>Silver Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Urgo Clean Ag®</b> <b>Manufacturer:</b> Urgo Ltd	<ul style="list-style-type: none"> <li>To be used under guidance from TVN, podiatry or vascular specialists.</li> <li>To reduce bioburden.</li> <li>Good for exuding wounds.</li> <li>Localised infection.</li> <li>Spreading infection.</li> <li>Systemic infection.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than 6 weeks at a time</li> <li>If wound is not showing signs of infection.</li> <li>Known sensitivity to silver.</li> <li>Patients undergoing MRI scan.</li> <li>Prior to radiotherapy.</li> <li>UrgoClean Ag is not suitable for use as a surgical sponge for heavily bleeding wounds.</li> <li>Do not use UrgoClean Ag in combination with hydrogen peroxide, organomercuric antiseptics or hexamidine.</li> <li>Where bioburden is not an issues.</li> </ul>	7 days	6cm x6cm 10cm x 10cm	10 £19.80 10 £4.51	<b>TVN/Pod only</b> <b>TVN/Pod only</b>
<b>UrgoTul Ag/Silver®</b> <b>Manufacturer:</b> Urgo Ltd	<ul style="list-style-type: none"> <li>To be used under guidance from TVN, podiatry or vascular specialists.</li> <li>To reduce bioburden.</li> <li>Localised infection.</li> <li>Spreading infection.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than 6 weeks at a time.</li> <li>If wound is not showing signs of infection.</li> <li>Known sensitivity to silver.</li> <li>Patients undergoing MRI scan.</li> <li>Prior to radiotherapy.</li> </ul>	7 days	10cm x12 cm	25 £81.00	<b>TVN/Pod only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

		<ul style="list-style-type: none"> <li>Where bioburden is not an issues.</li> </ul>				
<b>Acticoat Flex 7®</b> <b>Manufacturer:</b> Smith and Nephew	<ul style="list-style-type: none"> <li>To be used under guidance from TVN, podiatry or vascular specialists.</li> <li>To reduce bioburden.</li> <li>Localised infection.</li> <li>Spreading infection.</li> <li>Systemic infection.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time. Unless discussion with TVN</li> <li>If wound is not showing signs of infection.</li> <li>Known sensitivity to silver.</li> <li>Patients undergoing MRI scan.</li> <li>Prior to radiotherapy.</li> <li>Where bioburden is not an issue</li> </ul>	Up to 7 days	5cm x 5cm 10cm x 12.5cm	5 £32.40 5 £96.55	<b>TVN/Pod only</b> <b>TVN/Pod only</b>
<b>Aquacel Ag + Extra®</b> <b>Manufacturer:</b> ConvaTec	<ul style="list-style-type: none"> <li>To be used under guidance from TVN, not as first line antimicrobial.</li> <li>Moderately to highly exuding wounds.</li> <li>Ribbon to pack cavity wounds.</li> <li>To reduce bioburden.</li> <li>Localised infection.</li> <li>Spreading infection</li> <li>Systemic infection</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time.</li> <li>Dry wounds with minimal exudate.</li> <li>Known sensitivity to silver.</li> <li>Patients undergoing MRI scan.</li> <li>Prior to radiotherapy treatment.</li> <li>Where bioburden is not an issue.</li> </ul>	Up to 7 days	5cm x 5cm 10cm x 10cm	10 £20.20 10 £48.20	<b>PN, DN, NH, TVN/Pod</b> <b>PN, DN, NH, TVN/Pod</b>
<b>Chlorhexidine Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Bactigras® Manufacturer: Smith &amp; Nephew</b>	<ul style="list-style-type: none"> <li>Minor burns and scalds</li> <li>Lacerations, abrasions and other skin loss wounds</li> <li>Donor and recipient graft sites</li> <li>Where the wound is no more than 10% of the body area.</li> <li>Where there is a risk of infection, or on already infected wounds in conjunction with systemic antibiotics.</li> </ul>	<ul style="list-style-type: none"> <li>It should not be used in known cases of sensitivity or allergy to chlorhexidine.</li> <li>Flammable product – do not use near naked flames</li> </ul>	7 days	5cm x 5cm	50 £12.60	<b>TVN/POD only</b>
<b>DACC Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Cutimed Sorbact Swabs® Manufacturer: BSN</b>	<ul style="list-style-type: none"> <li>Dressing for contaminated, colonised or infected wounds.</li> <li>For deep, superficial, traumatic, postoperative or dehisced wounds.</li> <li>For ulcers (venous, arterial, diabetic, pressure).</li> <li>Suitable for fungal infections in the groin, skin folds or between digits.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use in combination with ointments and creams as the binding effect is impaired</li> </ul>	N/A	4cm x 6cm 7cm x 9cm	5 £8.70 5 £14.45	<b>PN, DN, TVN/Pod PN, DN, TVN/Pod</b>
<b>Iodine Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Inadine® Manufacturer: KCI/Acility</b>	<ul style="list-style-type: none"> <li>To be used on superficial infected wounds e.g traumatic wounds, superficial burns.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time.</li> <li>Use with caution in people with thyroid disorders.</li> <li>Not to be used in moderate to wet wounds.</li> </ul>	3 days	5cm x 5cm 9.5cm x 9.5cm	25 £8.25 10 £4.90	<b>DN, PN, NH, TVN DN, PN, NH, TVN</b>
<b>Iodoflex Paste® Manufacturer: S&amp;N</b>	<ul style="list-style-type: none"> <li>For treatment of chronic exuding wounds.</li> <li>For infected wounds.</li> </ul>	<ul style="list-style-type: none"> <li>Not to be used for more than two weeks at a time.</li> <li>To be used with caution in people who have thyroid disorders, pregnant or lactating women.</li> </ul>	7 days	5g 10g	5 £21.70 3 £25.98	<b>DN, PN, TVN DN, PN, TVN</b>
<b>Octenilin Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).

<b>Octenilin®</b> <b>Manufacturer:</b> Schulke	<ul style="list-style-type: none"> <li>Suitable for cleansing, moistening and decontamination of chronic skin wounds and to support the natural healing process.</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on individuals who have known allergic reactions to any of the contents. Not to be used on hyaline cartilage.</li> </ul>	3-5 minutes irrigation	350ml 20ml wound gel	1 £4.60 1 £4.78	<b>DN, PN TVN/POD only</b>
<b>Contact layer</b>						
<b>Dressings</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>UrgoStart Contact®</b> <b>Manufacturer:</b> Urgo Ltd	<ul style="list-style-type: none"> <li>Ulcers (low to moderately exuding leg ulcers, diabetic foot ulcers and pressure ulcers)</li> <li>Long standing acute wounds</li> </ul>	<ul style="list-style-type: none"> <li>Do not use on infected or critically colonised wounds</li> <li>Cancerous wounds</li> <li>Fistulas, which may reveal a deep abscess</li> <li>Known sensitivity to the dressing.</li> </ul>	7 days	5cm x 7cm	10 £30.70	<b>TVN/Podiatry only</b>
<b>N-A®</b> <b>Manufacturer:</b> 3M	<ul style="list-style-type: none"> <li>Ulcers (leg and pressure)</li> <li>Burns</li> <li>Other granulating wounds</li> <li>Allows free passage of exudate</li> <li>Requires secondary dressing</li> </ul>	<ul style="list-style-type: none"> <li>Do not use if allergic to silicone.</li> </ul>	7 days	9.5cm x 9.5cm	40 £14.00	<b>DN, PN, NH, TVN/Pod</b>
<b>Other podiatry dressings</b>						
<b>Other</b>	<b>Indications</b>	<b>Contraindications</b>	<b>Maximum wear time</b>	<b>Size</b>	<b>Cost per unit/pack size</b>	<b>Available to</b>
<b>Jelonet®</b> <b>Manufacturer:</b> Smith & Nephew	<ul style="list-style-type: none"> <li>Minor burns and scalds</li> <li>Donor and recipient graft sites</li> <li>Skin loss wounds</li> <li>Lacerations,</li> <li>Abrasions</li> <li>Leg ulcers</li> </ul>	<ul style="list-style-type: none"> <li>Flammable product – do not use near naked flames</li> </ul>	7 days	5cm x 5cm	50 £21.50	<b>TVN/Pod only</b>
<b>Cosmopore®</b> <b>Manufacturer:</b> Hartmann	<ul style="list-style-type: none"> <li>Sterile treatment of minor superficial wounds/injuries.</li> </ul>	<ul style="list-style-type: none"> <li>None</li> </ul>	2-7 days	5cm x 7.2cm	50 £2.70	<b>TVN/Pod only</b>

PN: Practice Nurse/DN: District Nurse/NH: Nursing Home/RH: Residential Home/TVN: Tissue Viability Nurse/Pod: Podiatry. All information and Drug Tariff prices correct as of May 21 This document was produced in conjunction with Frimley Health Foundation Trust, Frimley Health Community services and Virgin Care (now HCRG Care Group).